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Electric breast pump
Instructions for use

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ESPAÑOL
46-67

1. Important safety information

**⚠ Please read all instructions before using this product.
Keep these instructions for future reference.**

Warnings identify all instructions that are important for safety. Failure to observe these instructions can lead to personal injury or damage to the product. The following symbols and signal words show the significance of the warnings:

⚠ WARNING

Can lead to serious injury or death.

⚠ CAUTION

Can lead to minor injury.

NOTICE

Can lead to material damage
(not related to personal injury).

i Information

Useful or important information that is not related to safety.

When using electrical products, especially when children are present, basic safety precautions must always be followed. Please store the product where children can not reach.

⚠ WARNING

To avoid fire, electric shock, or serious burns:

- The breast pump should never be left unattended when plugged into a power source.
- Never operate an electrical device if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. If damage is found, immediately discontinue use and call Medela Customer Service.
- Do not use an electrical product that has been exposed to water or other liquids including:
 - Do not use while bathing or showering.
 - Never place or drop into water or other liquids.

– Do not run water over the breast pump.

- If an electrical product has been exposed to water or other liquids, do not touch and if possible unplug the device from electrical outlet, turn off and contact Medela.
- Do not expose the breast pump or its battery to excessive temperatures.
- Be aware of the risk of USB-C-socket of the pump or battery being shortcircuited by metal objects or fluids.
- Use only an external supply with the following specifications: DC 5V 2.5A. It must only be supplied with low voltage corresponding to the marking on the appliance.

⚠ WARNING

To avoid health risks and reduce the risk of severe injury:

- This is a single-user product. Use by more than one person may present a health risk and voids the warranty.
- Do not use the breast pump while operating a moving vehicle.
- Pumping may induce contractions or premature labour. Do not pump until after giving birth. If you become pregnant while breastfeeding or breast pumping, consult with a licensed healthcare professional before continuing.
- If infected with Hepatitis B, Hepatitis C, or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of

transmitting the virus to your baby through your breast milk.

- Wash all parts that come into contact with your breast and breast milk after every use.
- The tubing poses a risk of strangulation when being played with and small parts may be swallowed by toddlers.
- Close supervision is necessary when the breast pump or accessories are used in the vicinity of children.
- Keep this medical device away from magnetically susceptible medical devices such as cochlear implants, neurostimulators, stents and shunts.

⚠ CAUTION

To avoid fire or burns:

- Make sure the voltage of the power adaptor is compatible with the power source. See **chapter 17** for technical specifications.
- Always unplug electrical products immediately after charging.

- Do not place or store the product where it can fall into water.
- The breast pump and accessories are not heat-resistant. Keep away from heated surfaces or open flames.
- Do not use near flammable materials.

⚠ CAUTION

To avoid health risks and reduce the risk of injury:

- Before each use visually inspect the individual components for cracks, chips, tears, discoloration or deterioration. In the event that damage to the device is observed, please discontinue use until the parts have been replaced.
- This product cannot be serviced or repaired. Do not repair yourself. Do not modify the device.
- Never use a damaged device. Replace defective or worn parts.
- Use the breast pump only for its intended purpose/use as described in these instructions for use.
- Do not use the breast pump while sleeping or overly drowsy.
- If tubing becomes mouldy, discontinue use and replace tubing.
- Do not microwave or boil breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving (microwaving or boiling can also change the composition

of breast milk).

- Clean and disinfect all parts that come into contact with your breast and breast milk prior to first use.
- Only use Medela recommended parts with your breast pump.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain. If you are unsure about correct breast shield sizing please visit www.medela.com/fittingguide or see a lactation consultant/breastfeeding specialist who can help you get a proper fit.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful. See **chapter 7.1** for more information.

NOTICE

Take appropriate care in handling components:

- Collection cups and parts may become damaged if mishandled (e.g., dropped, or knocked over).

- Do not use the breast milk from bottles or components that show signs of damage.

2. Intended purpose

2.1 Indications for use

The Freestyle™ Mini Hands-free breast pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts. The Freestyle™ Mini Hands-free breast pump is intended for a single user. The breast pump is intended to be used in a home environment.

2.2 Contraindications

There are no known contraindications for the Freestyle™ Mini Hands-free breast pump.

3. Product description

3.1 Description of the device


The Freestyle™ Mini Hands-free breast pump is a personal-use electric breast pump that includes 2-Phase Expression technology and is capable of single and double pumping.

Mode of operation: Continuous. The operating life of the breast pump is 300 hours.

The expected operating life of the washable parts (breast shield, outer shell, yellow membrane, O-ring) is six months. The expected operating life of the washable parts (insert) is three months.

3.2 Your breast pump system features

Hands-free Collection Cups

1x Tubing

Short ends of tubing (connect to collection cups)

Long end of tubing (connects to motor unit)



Content may vary from picture.

4. Getting started

4.1 First steps

To be ready for expressing milk do the following before using your breast pump for the first time:


1. Remove the clip from the lanyard. Insert the clip's rope through the hole on the breast pump then push the clip into the loop. Reattach the clip to the lanyard.

2. **Charge the battery** of your Freestyle™ Mini Hands-free breast pump for three hours – see **chapter 4.3**.
3. **Separate all parts** that are exposed to the breast or breast milk prior to cleaning – see **chapter 5**.
4. Please ensure to boil all parts of the collection cups before first use – see **chapter 6.2** and **6.3**.

4.2 Powering your breast pump

By the rechargeable battery

Before first time use, charge the breast pump battery for three hours with the provided USB-C cable and power adaptor that meets the specifications.

→ The battery display shows all three bars when charging is complete.

For details on the battery display see **chapter 8.3**.

4.3 Charging the battery

⚠ WARNING

Only use the USB-C cable that comes with the breast pump. An inadequate power source or cable may result in a fire hazard, electric shock or malfunction equipment due to overvoltage. The warranty is only valid when using the Medela USB-C cable to charge your pump.

Charge the battery for three hours before using the breast pump for the first time.

1. Connect the USB-C cable to the motor unit.
 2. Connect the USB-C cable to the power adaptor (not provided).
 3. Plug the power adaptor (not provided) into a wall outlet.
- The battery display shows all three bars when charging is complete.

For details on the battery display see **chapter 8.3**.

⚠ WARNING

Do not expose the breast pump or battery to excessive temperatures.

- The breast pump cannot be used during charging
- Never charge the breast pump with a damaged USB power cable.
- Stop using the device immediately if you note any smoke or burning from the pump unit or power adapter.
- The breast pump cannot be turned on or operated when connected to the charger.

i Information

- The pump will not turn on when charging.
- The permissible ambient temperature for charging the battery is 5 °C to 40 °C (41 °F to 104 °F).


⚠ CAUTION

- Please fully charge the product before first use. Do not charge the battery pack for more than 4 hours. Dispose the waste battery packs under local laws and regulations.
- If the product is not used for a long time, the battery will run out of power. It needs to be plugged in for reactivation. After plugging in the power, wait for ten minutes for the screen to light up to indicate that it has been activated. Then continue charge the product till fully charged. If it is not activated for more than half an hour, please consult customer service.


5. Disassembling your collection cups

1. Open the collection cup using the tabs at the bottom.

2. Remove the yellow membrane from the breast shield.

6. Cleaning

It is important you do the following before using your breast pump for the first time: Disassemble all parts & clean by following your breast pump cleaning instructions.

Supplies needed:

- Dish washing liquid
- Clean soft brush
- Clean wash basin or dishwasher
- Drinking-quality water
- Clean pot for boiling water

Parts to wash or disinfect:

- Breast shields with O-ring
- Outer shells
- Yellow membranes
- Inserts

CAUTION


Separate and wash all parts that are exposed to breast milk immediately after use. This will help to remove breast milk residue and prevent growth of bacteria.

- Wash hands thoroughly with soap and water.
- Only use drinking-quality tap water or bottled water for cleaning.
- Do not place collection cup parts directly in the sink to rinse or wash. Use a clean wash basin used only for infant items.
- Do not use a dish towel to rub or pat items dry.
- Do not use disinfectants for cleaning.

NOTICE

Medela Freestyle™ Mini Hands-free breast pump has a closed system so milk cannot get into the tube when the pump is used as instructed. **Do not clean or rinse the tubing under any circumstance.**

- **Inspect the tubing before use. If you find condensation, breast milk and mould in the tubing, discontinue use and replace tubing.**
- Do not run the breast pump with wet tubing; doing so will cause damage to your breast pump.

NOTICE

Pay attention to the following:

- Take care not to damage parts of the collection cups during cleaning.
- Store the dry collection cups in a clean bag/container until next use.
- Do not store wet or damp parts as mould may develop.
- If you notice a white residue on your collection cups parts after boiling,

your water may have a high mineral content. Remove residue by wiping parts with a clean towel and allow to air dry.


- Distilled water is recommended when boiling parts to prevent substantial mineral build-up over time, which may compromise your parts.

6.1 Cleaning overview

Washing and disinfecting are two different activities. They must be done separately to protect you and your baby, and to maintain the performance of your breast pump.


Wash – To clean the surfaces of the parts by physically removing milk residues and other possible contaminants.

Disinfect – To efficiently lower the number of germs that may be present on the surfaces of the washed parts.

	Breast shield	Insert	Outer shell	Yellow membrane	O-ring
When to clean					
Before first use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Disinfect the parts according to chapter 6.3				
After each use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2				
Once per day after use	Disassemble the parts according to chapter 5 Wash the parts according to chapter 6.2 Disinfect the parts according to chapter 6.3				
When to clean	Tubing 		Motor unit 		
As needed	Do not clean or rinse the tubing under any circumstance.		Wipe with clean, damp towel. Refer to chapter 6.4.		

6.2 Washing

Rinse and wash by hand (before first use and after each use)


Do not place the parts directly in the kitchen sink for rinsing and washing. Use a dedicated wash basin for infant items.

1. Rinse the disassembled parts, except for the tubing, with cold, clear drinking quality water (approx. 68 °F/20 °C).
2. Clean these parts with plenty of warm, soapy water (approx. 86 °F/30 °C). Use a commercially available dish soap and clean soft brush.
3. Rinse the parts with cold, clear drinking quality water for 10 to 15 seconds (approx. 68 °F/20 °C).
4. Allow to completely dry after washing.

or

Wash in dishwasher (maximum once per day)


1. Place the disassembled parts, except for the tubing, on the top rack or in the cutlery section. Use a commercially available dishwashing detergent.
2. Rinse the yellow membrane under running water after dishwashing.
3. Allow to completely dry after washing.

Information

If using the dishwasher, parts may become discolored. This will not impact the functioning of the parts.
Do not wash the tubing.


Drying

1. Allow to air dry on a clean, unused dish towel or paper towel

6.3 Disinfecting

Disinfect on stovetop (before first use and once per day after use)

1. Cover the disassembled, washed parts, except for the tubing, with plenty of water and boil at least for 10 minutes. Ensure that the parts are not in direct contact to the bottom of the pot.
2. Allow to dry after disinfecting.

or

Disinfect in microwave (once per day after use)

1. Use Quick Clean™ Micro-Steam™* bags in the microwave in accordance with the instructions on the bags.
2. Allow to dry after disinfecting.

* Refer to local website/shops for availability in your country.

Drying

1. Allow to air dry on a clean, unused dish towel or paper towel.
 2. Put the clean and dry parts in a clean storage bag or a clean environment.
- Do not store parts in an airtight container/bag if moist. It is important that all residual moisture dries.

6.4 Cleaning the motor unit


Supplies needed:

- Drinking-quality water
- Dish washing liquid
- Clean dish towel or paper towel

Parts needed:

- Motor unit

CAUTION

To prevent electric shock and damage to the material:

- Unplug the breast pump before cleaning the motor unit.
- Do not immerse the motor unit in water; do not run water over the motor unit.

NOTICE

- Do not spray or pour liquid directly onto the pump.
- When cleaning the motor unit display, only use a damp, soft towel with mild washing-up liquid. Other cleaners or abrasive cloths may degrade the display performance.





1. Switch the breast pump off.
2. Unplug the breast pump from the power source.
3. Wipe the motor unit with a clean towel moistened with drinking-quality water.
4. Wipe the motor unit dry with a clean, dry towel

7. Assembling your breast pump

Only use clean and dry parts for assembly. Refer to **chapter 6** for cleaning.

7.1 Selecting the correct size of the breast shield and insert

Make sure to use the correct breast shield size. It is an essential component of effective pumping, helping to optimise milk flow.

Your nipple diameter (mm)				
11 mm or smaller	12 mm, 13 mm or 14 mm	15 mm, 16 mm or 17 mm	18 mm, 19 mm or 20 mm	21 mm or bigger
				
				
Insert 15 mm	Insert 18 mm	Insert 21 mm	Breast shield 24 mm	
				
Included		Breast shield 21 mm	Breast shield 24 mm	Breast shield 27 mm
				
Accessories to be purchased				

1. Select the breast shield size by measuring the diameter of your nipple.

2. Select the breast shield size that is slightly larger than your nipple.

2 x 24mm size breast shields and 2 x 3 inserts of 15, 18 and 21mm are included in the pack.

For more information on breast shield sizing and other available sizes, visit www.medela.com.

7.2 Assembling the collection cups

Only use clean and dry parts for assembly. Refer to chapter 6 for cleaning.

Parts needed:

- Tubing
- Breast shield with O-ring
- Motor unit
- Yellow membrane
- Outer shell
- Inserts

CAUTION

To prevent contamination of your milk:

- Wash hands thoroughly with soap and water before touching the motor unit, parts of the collection cups and breasts.
- Dry your hands with a fresh towel or a single-use paper towel.
- Avoid touching the inside of the collection cups components.


To make sure that your breast pump works properly and safely:

- Check collection cups components for wear or damage before use. Replace if necessary.
- Always inspect all parts prior to use for cleanliness. If dirty, clean according to **chapter 6**.
- Only use genuine Medela parts.


NOTICE

To prevent damage to the breast pump all components must be completely dry before use.


10


1. Wash hands thoroughly with soap and water before touching the motor unit, parts of the collection cups and breasts. Dry your hands with a clean towel.


2. Confirm that the O-ring is applied to the breast shield.

3. First assemble the valve of the yellow membrane on the front of the breast shield.


4. Align the disc of the yellow membrane on the designated circular slot using the alignment tab at the top to ensure that it is at the right place. Press around the edges to make sure it is fully assembled.


5. Insert the pouring spout into the opening on the shield.


6. Push the outer shell and breast shield together tightly, ensuring it is fully sealed all the way around until you hear a snapping sound.

7. If needed add the insert according to the sizing table.

8. Connect the long end of the tubing to the motor unit as far as it will go.

9. Connect the short end of the tubing into the tubing connection port on the collection cup.


i Information

O-ring fits all breast shield sizes. To disassemble and assemble the O-ring to a different size of breast shield follow the guidance below.

- **To disassemble the O-ring:** Start by pulling out the bottom tab (as shown in the left image below) from the breast shield.
- **To reassemble the O-ring:** Start from the top (as shown on the right image below) of the breast shield, stretch it over the outer edge of the breast shield. Ensure it goes on with the ridges facing outwards. Press along the groove, checking that it is in place all the way around.

To disassemble:


Pull out the bottom tab from here

To assemble:


Start placing the O-ring from here.

Ensure that the
ridges are
straight

7.3 Preparing for single pumping

1. Assemble one collection cup as described in **chapter 7.2**.

Insert the unused tubing end into the Y-shaped tubing holder.
Correctly assembled system (for single pumping).

2. Place the collection cup in your bra according to **chapter 7.5**. Start pumping as described in **chapter 9**.


7.4 Preparing for double pumping

i Information

- Double pumping saves time and increases the nutritional value of the milk.


1. Assemble the two collection cups as described in **chapter 7.2**.

Correctly assembled system (for double pumping).


2. Place the collection cups in your bra according to **chapter 7.5**. Start pumping as described in **chapter 9**.


7.5 Placing collection cups within your bra


1. Ensure that the nipple is centered in the tunnel of the breast shield.


2. The hands-free collection cups can be worn under your regular nursing bra. Your bra will hold the collection cup in place. If you find it uncomfortable or painful during use, try switching to a looser, more stretchy nursing bra.

3. The tubing can be connected facing inwards (towards the breasts) or outwards (towards the arms) based on personal preference. Ensure that the tubing is not pinched or bent during use.

4. The hands-free collection cups can be used in a vertical or slightly tilted position for better visibility of the nipple and milk ejection. The scale on the collection cup gives you an approximate indication of the volume during pumping.

5. With the lanyard the breast pump can be worn cross-body or on the waist.


8. Operating your breast pump

8.1 Terms and designations


2-Phase Expression	 Through research, Medela found that babies nurse in 2 phases – Stimulation and Expression. This research is the basis for the technology in all Medela breast pumps.
Stimulation phase	Fast sucking/pumping mode to stimulate milk flow.
Expression phase	Slower sucking/pumping mode for gentle and efficient milk removal after milk has started flowing.
Maximum comfort vacuum	Highest vacuum setting where pumping still feels comfortable. Different for every mother.
Let-down	The reflex that makes the milk begin to flow from the breast.

8.2 Operating and display elements

Control elements

Display elements

8.3 Power and battery status display

	The battery display shows the state of the battery and its charging level. → The display of all three bars indicates that the battery is full.
	With decreasing battery capacity fewer bars appear. → Two bars: enough for approx. 90min operation. → One bar: enough for approx. 70 min operation.
 (flashing)	→ When the breast pump's last single bar is flashing, it indicates that the battery is almost empty. The pump has to be charged.

9. Expressing your breast milk

⚠ CAUTION


- Do not try to express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful.

14

NOTICE

Make sure that the tubing is not kinked or pinched while pumping.

1. Make sure that the collection cups are connected and positioned properly.

2. Long press the On/Off/Pause button  on the front until the lights on the screen turn on, then release the button.

→ The breast pump begins in the Stimulation phase.


3. Adjust the vacuum with the Increase vacuum () and Decrease vacuum () buttons to find a comfortable level in Stimulation phase.

4. As soon as your milk begins to flow, press the Let-down button  to go to the expression phase.

→ The breast pump proceeds with the Expression phase.

→ If the Let-down button  is not pressed within one minute, your breast pump will automatically go to Expression phase.

5. Adjust the vacuum with the Increase vacuum () and Decrease vacuum () buttons to find your Maximum Comfort Vacuum in Expression phase.


Finding your Maximum Comfort Vacuum

Maximum Comfort Vacuum is the highest vacuum level where pumping still feels comfortable.

1. Once you are pumping in the Expression Phase, increase the vacuum with the () button until pumping feels slightly uncomfortable (not painful).
2. Then decrease the vacuum slightly with one press of the () button.

i Information

- The pump is not designed to be used while lying down or sleeping. However, you can slightly lean forward or backward during your pumping session.
- The maximum volume/capacity for the collection cup is 150 ml/5 oz. When this limit is reached, transfer the milk to a bottle before continuing to express.

6. When your pumping session is over, press the On/Off/Pause button  to stop the pumping operation.
→ The breast pump stops automatically after 30 minutes if there is no user interaction with the pump. The breast pump will automatically switch off after 10 minutes of pausing.
7. Long press the On/Off/Pause button  on the front to switch the breast pump off.

Pausing

In both stimulation and expression phase the breast pump can be paused.

This allows you to readjust yourself or your collection cups, or attend to your surroundings.

1. For pausing, press the On/Off/Pause button  while the breast pump is running.
2. To resume pumping, press the On/Off/Pause button  again.

If you want to start a completely new pump session, switch the breast pump off. Long press On/Off/Pause button  to start again.

i Information


- This breast pump will automatically turn off after working 30 minutes. The breast pump will automatically switch off after 10 minutes of pausing.
- To resume pumping, you will need to start a new session.
- If you need to pause for more than two minutes, we recommend restarting with stimulation phase.

10. Handling of breast milk


10.1 Removing the collection cups from your bra

i Information


To avoid milk spilling out of the tunnel, when finished, keep the pump running and slightly lean forward to ensure all the milk inside the breast shield tunnel empties into the collection cups.

1. While the pump is still running, slightly lean forward and ensure all milk from the nipple tunnel enters the collection container. Turn off pump.


2. Detach the tubing from the collection cups.

3. Carefully remove the collection cups from your bra.

10.2 Pouring milk from collection cups

1. Place the collection cups in an **upright position** on a surface where they will not be easily knocked over. You can use the volume scale markings on the shield to read the collected milk volume. The scale on the outer shell gives you an approximate indication of the volume during pumping.

2. **With the rounded outer shell facing down**, pour breast milk through the pouring spout into a bottle or container for storage. **Do not use collection cups to store milk. Do not store assembled cups in the refrigerator.**

3. Fully disassemble and clean after each use and disinfect once per day as described in **chapter 5** and **chapter 6**.

10.3 General storage guidelines

General breast milk handling guidelines

- Write the date on the storage container. Include your child’s name if you are giving the milk to a childcare provider.
- Gently swirl the container to mix the cream part of the breast milk that may rise to the top back into the rest of the milk. Do not shake the milk. This can make some of the milk’s valuable parts break down.
- Refrigerate or chill milk right after it is pumped, if possible. You can put it in the refrigerator, place it in a cooler or insulated cooler pack, or freeze it in small (2- to 4-ounce) batches for later feedings. Pumped milk, at room temperature ($\leq 77^{\circ}\text{F}$) is safe to use without being refrigerated for up to 4 hours after pumping.

Human milk storage guidelines

Storage locations and temperatures			
Type of breast milk	Countertop 77°F (25°C) or colder (room temperature)	Refrigerator 40°F (4°C)	Freezer (-18°C) or colder
Freshly expressed or pumped	Up to 4 Hours	Up to 4 Days	Within 6 months is best Up to 12 months is acceptable
Thawed, previously frozen	1–2 Hours	Up to 1 Day (24 hours)	Never refreeze human milk after it has been thawed
Leftover from a feeding (baby did not finish the bottle)	Use within 2 hours after the baby is finished feeding		


*Recommended storage times are important to follow for best quality.

For information on storing your breast milk, visit <https://www.medela.com/en/breastfeeding-pumping/articles/pumping-tips/how-to-store-freeze-and-thaw-breast-milk>.

Learn more about storing your breast milk at <https://www.cdc.gov/breastfeeding/breast-milk-preparation-and-storage/handling-breastmilk.html> and <https://womenshealth.gov/breastfeeding/pumping-and-storing-breastmilk>.

11. Maintenance and care

11.1 Power adaptor care

NOTICE

To keep your power adaptor and USB-C cable in good working condition, pay attention to the following:
Do not unplug the cable or the power adaptor by pulling on the cord, but hold it by the connector housing and adaptor body.

Traveling internationally

i Information

While traveling, use the USB-C cable provided and pair with the appropriate USB power source for your current location.

11.2 Battery maintenance

Your breast pump contains a built-in (not replaceable) lithium-ion rechargeable battery.

CAUTION

- Please use the charging USB-C cable provided with this product. Do not plug the USB-C cable of any other manufacturer into the motor for use, so as not to cause damage to the motor.
- Please use adaptor, select an adaptor with input voltage AC 100–240V and output DC 5V and 2.5 A.

NOTICE

Be aware of the risk of USB-C-socket of the pump or battery beeing shortcircuited by metal objects or fluids.

Traveling internationally


Information

This breast pump has a lithium-ion battery. Carrying such batteries on airplanes may be restricted by the country you are visiting. Please consult with the country you are visiting to find out if there are any restrictions that pertain to traveling with lithium-ion batteries. Be aware that the performance of (product name) may be affected by the external environment, in view of the uncertainty and instability of the travel environment.

11.3 Long-term storage

For long term storage keep the pump in a dry and cool place.


Before storing your breast pump for an extended amount of time remember to:

☐ Clean the collection cup
(refer to **chapter 6**)

☐ Clean the motor unit
(refer to **chapter 6.4**)

☐ Make sure that the battery
is charged (see below)

It is best to fully charge the battery before storing your breast pump for an extended period of time. Refer to **chapter 8.3** for information on the battery charging state.

CAUTION

- If the product is not used for a long time, the battery will run out of power. It needs to be plugged in for reactivation. After plugging in the power, wait for ten minutes for the screen to light up to indicate that it has been activated.

12. Troubleshooting

In case of an unexpected behaviour of your breast pump check with the troubleshooting table if you find the issue in the column "Problem" and follow the instructions in the column "Solution".

Troubleshooting table

Problem	Solution
The breast pump does not enter the standby mode after you pressed the On/Off/Pause button	- You pressed the On/Off/Pause button for too short or too long a time. - Press the On/Off/Pause button until you hear an audio signal, then release it. If the pump is in storage mode, it will take somewhat longer (approx. one second) until the audio signal is emitted.
The display does not light up after you pressed the On/Off/Pause button	- If your pump has not been used in a while and it does not start as usual, it might be in hibernation mode. - To reset your breast pump, press and hold the On/Off/Pause button for seven to ten seconds. - If there is no change, make sure that the breast pump is connected to a properly working power supply. - If there is still no change, contact the Medela Customer Service.
The timer does not count up after you pressed the On/Off/Pause button	Contact the Medela Customer Service.
The breast pump generates no vacuum (motor not working) after you pressed the On/Off/Pause button	- Make sure that the breast pump is switched on. - Contact the Medela Customer Service.
The breast pump generates only constant vacuum instead of sucking cycles	Contact the Medela Customer Service.
The breast pump stops operation or cannot be turned on (emits an audio signal and all elements of the battery display are flashing)	The breast pump refuses operation due to battery over-heating. Allow the motor unit to cool down, then restart the breast pump.
All elements of the battery display are flashing when you are trying to charge your breast pump	Charging is not possible when the battery temperature is too high. Allow the motor unit to cool down, then resume charging.
There is low or no suction	- Make sure that all collection cups components are clean and dry and that connections are secure. - While pumping, make sure the breast shields form a complete seal around the breast. - When single pumping, make sure that the unused tubing end is correctly plugged into the tubing holder. - If suction does not improve after having followed these steps, contact Medela Customer Service.
The breast pump is not responding as expected	- To reset your breast pump, press and hold the On/Off/Pause button for seven seconds. - If there is no change, contact the Medela Customer Service.
The breast pump exterior got wet	- Unplug the breast pump from the power source and switch it off. - Dry off the exterior of the breast pump.
The breast pump has been submerged in water	- Unplug the breast pump from the power source and switch it off. - Contact Medela Customer Service.

If you have not resolved the problem with your breast pump or you have further questions, please contact Medela Customer Service. For contact data visit www.medela.com. Under "Country" choose your country.

13. Warranty

This product is warranted by Medela LLC to the original retail purchaser to be free from defects in material and workmanship for the period of 1 year for pump mechanism and 90 days for other components, from the date of purchase. Warranty can only be claimed in the country of purchase. In the event of a covered defect, Medela will, at its option, replace or repair the defective component, without charge. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO A DURATION OF 1 YEAR FROM DATE OF PURCHASE. SOME AREAS DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE OR IN YOUR COUNTRY.

Before you make a claim under this warranty, it may save time and expense to contact Medela Customer Service (see contacts below). You may also contact us for additional information concerning this warranty. If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, together with your proof of purchase and a brief statement of the problem to the following address:

Medela LLC – Returns, Door 4501
1101 Corporate Dr., McHenry, IL 60050, USA
ATTENTION: RETURNS
Medela Customer Service (toll free): 1-800-435-8316
customer.service@medela.com
<https://www.medela.com/en-us/breastfeeding-pumping/services/customer-service>

Contact Customer Service for a return authorization number and a prepaid shipping label before returning your product. Returns are not accepted without an authorization number. Find information on the international warranty on www.medela.com/ewarranty.

14. Disposal


Disposal of your breast pump

At the end of its operating life, separate the parts of your breast pump and dispose of according to the following instructions.

Breast shields, outer shells, tubing, O-rings, membranes and inserts

The parts are made of plastics that are not harmful to the environment when disposed of as household waste. Recycle or dispose of according to local regulations.

Motor unit and power adaptor

Do not dispose of electric or electronic equipment together with unsorted municipal waste, but collect it separately.

In the European Union the manufacturer or its vendor must take back waste equipment. Other countries may have similar collection and recycling systems. Dispose of the waste equipment according to local regulations.

Inquire at the point of sale or contact your local authority for appropriate collection points for waste equipment. The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment.

CAUTION

Hazardous materials present in electric/electronic equipment may harm human health and the environment if the waste equipment is not disposed of properly.

15. Meaning of symbols

The following tables explain the meaning of the symbols found on the product parts and its packaging.

	General safety alert symbol, points to information related to safety. ¹		Defines a temperature range [e.g. for operation, transport or storage]. ¹⁰
	Defines a relative humidity range (e.g. for operation, transport or storage). ²		Defines an atmospheric pressure range (e.g. for operation, transport or storage). ¹¹
	Indicates the location of the On/Off button. ³		Indicates the manufacturer. ¹²
	Indicates the date of manufacturing. ⁴		Indicates the part number of the device. ¹³
	Indicates the serial number of the device. ⁵		Read and follow the instructions for use. ¹⁴
	Indicates compliance with international requirements for protection from electric shock (Type BF applied parts). ⁶		Do not dispose of electric/electronic devices together with unsorted municipal waste (dispose of the device in accordance with local regulations). ^{15*}
	The CE mark indicates conformity with the applicable European directives.*		Indicates the degree of protection against ingress of foreign objects and moisture. ¹⁶
	Contains fragile goods. Handle with care. ⁷		Keep away from sunlight. ¹⁷
	Keep away from rain. Keep in dry conditions. ⁸		Indicates the compliance with the requirements of the Federal Communications Commission. ¹⁸
	Indicates the lot/batch number. ⁹		Indicates that the package is capable of being recycled. ¹⁹
	The packaging contains products intended to come in contact with food according to Regulation (EC) No 1935/2004.		Indicates the material from which an item is made. ^{20*}
	Indicates the model number of the device. ²¹		

References

- 1 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.2 Symbol 2 General warning sign.
- 2 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.8 Humidity Limit / ISO 7000-2620, Graphical symbols for use on equipment, Humidity Limit.
- 3 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 29 Stand-by.
- 4 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.3 Manufacturing Date / ISO 7000-2497, Graphical symbols for use on equipment, Date of manufacture.
- 5 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.7 Serial number / ISO 7000-2498, Graphical symbols for use on equipment, Serial number.
- 6 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 20 Type BF applied parts.
- 7 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.1, Fragile, handle with care / ISO 7000-0621, Graphical symbols for use on equipment, Fragile, handle with care.
- 8 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.4, Keep away from rain / ISO 7000-0626, Graphical symbols for use on equipment, Keep away from rain.
- 9 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.5 Batch code / ISO 7000-2492, Graphical symbols for use on equipment, Batch code.
- 10 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.7 Temperature Limit / ISO 7000-0632, Graphical symbols for use on equipment, Temperature Limit.
- 11 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.9 Pressure Limit / ISO 7000-2621, Graphical symbols for use on equipment, Pressure Limit.
- 12 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.1 Manufacturer.
- 13 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.6 Article number / ISO 7000-2493, Graphical symbols for use on equipment, Catalogue number.
- 14 IEC 60601-1, Medical electrical equipment – Part 1: General requirements for basic safety and essential performance, Table D.2 Symbol 10 Refer to instruction manual/booklet.
- 15* EN 50419, Marking of Electrical and Electronic Equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE).
- 16 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.3 Symbol 2 IP Code IEC 60529: Degrees of protection provided by enclosures (IP Code).
- 17 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.2 Keep away from sunlight / ISO 7000-0624, Graphical symbols for use on equipment, Keep away from sunlight.
- 18 Code of Federal Regulations, Title 47, Part 15b.
- 19 ISO 7000-1135, Graphical symbols for use on equipment, General symbol for recovery/recyclable.
- 20* 97/129/EC, Commission Decision establishing the identification system for packaging materials.
- 21 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.10 Model number

* These symbols are not applicable for the US, Canada, Mexico market.

16. International regulations

16.1 Electromagnetic compatibility (EMC)

The Freestyle™ Mini Hands-free breast pump is intended to be used in a home environment and is EMC-tested in conformity with the requirements of the standard IEC 60601-1-1-2.

The breast pump needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in these instructions for use. Portable and mobile RF communications can affect the breast pump.

 **WARNING**

To prevent malfunction resulting from electromagnetic interference:

- The electric breast pump Freestyle™ Mini Hands-free should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the electric breast pump Freestyle™ Mini Hands-free should be observed to verify normal operation in the configuration in which it will be used.
- Portable wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies, RFID, other electrical equipment and also including peripherals such as antenna cables and external antennas can affect the electric breast pump and should be kept at a distance of at least 30 cm away from the breast pump including cables. Otherwise, degradation of the performance of this equipment could result.
- Use of accessories or cables other than those provided by Medela of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

NOTICE

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Guidance and manufacturer's declaration – electromagnetic emissions

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Emission tests	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	The breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment. The Freestyle™ Mini Hands-free breast pump has no essential performance but was tested for immunity to electromagnetic disturbances and passed using the following criteria:

1. No change in mode and operation of the breast pump.
2. The pump can be controlled by changing the vacuum level and/or pressing the On/Off button.

Reset (The following sequence can be used to reset the device in the event of an EM Disturbance)

1. Unplug the USB-C cable from the motor unit.
2. Unplug the power adaptor from the wall outlet.
3. Attach a breast pump set to the tubing port (if none is attached).
4. Press and hold the On/Off button for seven to ten seconds.
5. Plug the power adaptor back into the wall outlet.
6. Plug the USB-C cable into the supply socket of the motor unit.
7. Press the On/Off button to switch the breast pump on.
8. Press the Start/Pause button to begin pumping.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact discharge +/- 15 kV air discharge	+/- 8 kV contact discharge +/- 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 5%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV, ± 1 kV Line-to-line	± 0,5 kV, ± 1 kV Line-to-line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U _T for 1 cycle at 0° 70 % U _T for 25 cycles at 50 Hz at 0° for 30 cycles at 60 Hz at 0° 0 % U _T for 250 cycles at 50 Hz at 0° for 300 cycles at 60 Hz at 0°	0 % U _T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0 % U _T for 1 cycle at 0° 70 % U _T for 25 cycles at 50 Hz at 0° for 30 cycles at 60 Hz at 0° 0 % U _T for 250 cycles at 50 Hz at 0° for 300 cycles at 60 Hz at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the breast pump requires continued operation during power mains interruptions, it is recommended that the breast pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	It may be necessary to position the breast pump further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Proximity magnetic fields IEC 60601-1-2, Table 11 (IEC 61000-4-39)	30 kHz, 8A/m, CW modulation 134.2kHz, 65A/m, Pulse modulation 2.1kHz 13.56MHz, 7.5A/m, Pulse modulation 50kHz	30 kHz, 8A/m, CW modulation 134.2kHz, 65A/m, Pulse modulation 2.1kHz 13.56MHz, 7.5A/m, Pulse modulation 50kHz	Mains power quality should be that of a typical commercial or hospital environment.

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Table of frequencies of portable and mobile transmitters for which the recommended separation distance is 30 cm (12 inches):

Band (MHz)	Service
380 - 390	TETRA 400
430 - 470	GMRS 460, FRS 460
704 - 787	LTE Band 13, 17
800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5
1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS
2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7
5 100 - 5 800	WLAN 802.11 a/n

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Freestyle™ Mini Hands-free breast pump including cables specified. Otherwise, degradation of the performance of this equipment could result.

17. Technical specifications


Summary of important technical specifications

Vacuum performance
Rated atmospheric Pressure: 850HPa-1060HPa
-50 to -270 mmHg, ± 20mmHg (-66 to -360 HPa, ±26 HPa)
45 to 111 cpm, ±5 cpm

Input	100 to 240 V~ 50/60 Hz max. 0.4 A
Output	5 V --- 2.5 A 12.5 W
No-load power consumption	Lower than 0.1 W

Motor unit


Battery capacity and type	2200mAh Li, Size 18650
Size	113 x 68 x 42 mm, ±1 mm
Weight	205 g, ±10 g

Operating temperature range
41°F to 104°F (5°C to 40°C)

Operating/Transport/Storage
humidity range
(20% to 80% relative humidity)

Ambient pressure range
(850HPa to 1060HPa)

Transport/storage
temperature range
-4°F to 131°F (-20°C to 55°C)

Materials touching skin or coming in contact with milk

- Breast shield : Polypropylene
- Outer shell : Polypropylene, Thermoplastic Elastomer
- Membrane : Silicone
- Insert: Silicone

Made without BPA.

Medela wordmark and logo and Freestyle are registered in the U.S. Patent and Trademark Office and the Canadian Intellectual Property Office. The Science of Care, Turning Science into Care and Medela Family are trademarks of Medela.

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Referencias

- 1 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.2 Symbol 2 General warning sign.
- 2 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.8 Humidity Limit / ISO 7000-2620, Graphical symbols for use on equipment, Humidity Limit.
- 3 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 29 Stand-by.
- 4 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.3 Manufacturing Date / ISO 7000-2497, Graphical symbols for use on equipment, Date of manufacture.
- 5 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.7 Serial number / ISO 7000-2498, Graphical symbols for use on equipment, Serial number.
- 6 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.1 Symbol 20 Type BF applied parts.
- 7 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.1, Fragile, handle with care / ISO 7000-0621, Graphical symbols for use on equipment, Fragile, handle with care.
- 8 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.4, Keep away from rain / ISO 7000-0626, Graphical symbols for use on equipment, Keep away from rain.
- 9 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.5 Batch code / ISO 7000-2492, Graphical symbols for use on equipment, Batch code.
- 10 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.7 Temperature Limit / ISO 7000-0632, Graphical symbols for use on equipment, Temperature Limit.
- 11 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.9 Pressure Limit / ISO 7000-2621, Graphical symbols for use on equipment, Pressure Limit.
- 12 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.1 Manufacturer.
- 13 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.6 Article number / ISO 7000-2493, Graphical symbols for use on equipment, Catalogue number.
- 14 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.2 Symbol 10 Refer to instruction manual/booklet.
- 15* EN 50419, Marking of Electrical and Electronic Equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE).
- 16 IEC 60601-1, Medical electrical equipment – Part 1: General Requirements for basic safety and essential performance, Table D.3 Symbol 2 IP Code IEC 60529: Degrees of protection provided by enclosures (IP Code).
- 17 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.3.2 Keep away from sunlight / ISO 7000-0624, Graphical symbols for use on equipment, Keep away from sunlight.
- 18 Code of Federal Regulations, Title 47, Part 15b.
- 19 ISO 7000-1135, Graphical symbols for use on equipment, General symbol for recovery/recyclable.
- 20* 97/129/EC, Commission Decision establishing the identification system for packaging materials.
- 21 ISO 15223-1, Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements, Clause 5.1.10 Model number

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