



IMPORTANT SAFETY INFORMATION: NOVA ORTHO-MED INC. BED RAILS RECALL

Background Information: The U.S. Consumer Product Safety Commission (CPSC) and NOVA Ortho-Med Inc. are announcing the voluntary recall of two models of the company's Bed Rails. This voluntary recall involves the Home Bed Rail (model 6093), and Home Bed Rail with Legs (model 6094).

Hazard: When attached to an adult's bed, users can become entrapped within the bed rail or between the bed rail and the side of the mattress. This poses a serious entrapment hazard and risk of death by asphyxiation.



NOVA Item 6093



NOVA Item 6094

Remedy: Consumers should immediately stop using the recalled bed rails and use extra caution when getting in and out of bed if a replacement bed rail is not used. For consumers who own the Home Bed Rail (model 6093), and/or Home Bed Rail with Legs (model 6094) contact NOVA Ortho-Med Inc. for a replacement bed rail unit for item 6093 and a repair kit for item 6094.

Please contact our customer service department at 1-800-557-6682 from 7:30a.m. to 5:30p.m. PST, Monday through Friday or email us info@novajoy.com



To verify if this recall applies to your product, please inspect your product for a sticker indicating model number “6093” or “6094” with the NOVA brand logo. If you are facing the product this sticker will be located on the left side of the handle. If you have a bag currently attached to your unit, you will need to remove the bag to locate the sticker.

For owners with Home Bed Rail (model 6093): You will be asked to verify the model number of your unit. You will then be asked to provide your address so we are able to send 2x sets of “RECALLED- DO NOT USE” sticker which you will apply to your product. After the sticker is applied, please contact our customer service number (800)-557-6682 to confirm you have applied the sticker defacing the product, then you will be instructed to discard the recalled unit into the trash and we will send you a replacement CSPC approved 6094 unit.

For owners with Home Bed Rail with Legs (model 6094): You will be asked to verify the model number of your unit. You will then be asked to provide your address so we are able to send a free repair kit.



NOVA Item 6093



NOVA Item 6094