

# Innovative Technology. Easy-to-Use. Clinically Effective.

The original Nisus pump is portable and delivers precise pressure for optimal power and performance to provide a universal solution for all wounds across the patient's full continuum of care.



The Cork Medical Nisus Negative Pressure Wound Therapy System is designed to provide gentle powered suction to treat and promote wound healing by aiding in the removal of excess exudates, infectious material, and tissue debris.

## NISUS NPWT SYSTEM FEATURES

The Cork Medical Nisus Negative Pressure Wound Therapy System is indicated for use in patients who would benefit from negative pressure wound therapy particularly as the device may promote wound healing by the removal of excess exudates, infectious material, and tissue debris.

- ▶ Recommended for patients with chronic, acute, traumatic, sub-acute and dehisced wounds, partial-thickness burns, pressure ulcers, diabetic ulcers, and venous ulcers, flaps, and grafts
- ▶ FDA approved for use in hospitals, long term care facilities, and home-care settings
- ▶ Light-weight and portable, with up to 24hrs of battery life (can extend to 48hrs with sleep mode)
- ▶ User friendly interactive interface with easy to read LCD screen
- ▶ Step-by-step troubleshooting guide for alarms such as Canister Full, Vacuum Leak, and Blockage
- ▶ Multiple unique features for clinical enhancement of negative pressure wound therapy such as: patient runtime, sleep feature, alarm log, lock restriction, Variable Intermittent and Continuous Therapy Modes



FDA Approved For Use In  
Hospitals, Long-Term Care Facilities,  
And Home-Care Settings



Sleep Feature Technology  
Allows Users To Dim Light For  
Long Lasting Battery Life



Multiple Features For  
Clinical Enhancements



Universal Solution For  
All Wound Types



Customizable Therapy Modes



Dynamic Troubleshooting Guide





## COMPONENT SPECS

Comprised of three essential components that actively work together to promote wound healing, the Nisus NPWT System includes the Nisus NPWT Pump, the Nisus NPWT Canisters, and the Cork NPWT Wound Dressing Kit. Please use under the advisement of a physician. Carefully read all therapy and device instructions and safety information prior to use.

Use only Cork Medical components with the Nisus Negative Pressure Wound Therapy System.

- ▶ Battery charger and carrying bag included

## UNIVERSAL NISUS PLATFORM & ACCESSORIES



**Nisus NPWT Pump**  
CMPP-100



**Nisus NPWT Canisters**  
CPC-250 (250ML)  
CPC-500 (500ML)




**NPWT Wound Dressing Kit**  
WDK-250 (SMALL)  
WDK-300 (MEDIUM)  
WDK-400 (LARGE)



**NPWT Cleanse Therapy Kit**  
CLEANSE-KIT-A-18



**Bridge Foam**  
NIS-BRDG-A-16



**Optifoam White Foam**  
FOAM-WHITE-A-16



**24" Extension Tubing**  
CSP-100-ES-EXT



**Y-Connector**  
CWY-100

