

medela 

Freestyle™ Mini Hands-free

Big on Performance.
Light on You.

Turning Science into Care



Discover Medela's ultra small and easy-to-carry, Freestyle™ Mini Hands-free breast pump. The essential hands-free pumping experience, recommended by 9 out of 10 moms.* The palm-sized motor and collection cups (only 2.7oz each) ensure a comfortable pumping experience without added weight, while the anatomically shaped cups reduce pressure on sensitive breasts and fit discreetly in your bra. Perfect for on-the-go moms.



Ultra small & light



2-Phase Expression technology
mimics babies' natural
sucking rhythm



Digital display
Intuitive control with
4 buttons and 9 pre-
programmed settings



Built-in battery
Lasts up to 2 hours on
a single charge



USB-charging
For powering up
anywhere you are



3 Parts**
Easy to clean and
assemble



Anatomically designed
Fits most breast shapes
and nursing bras



Capacity
of 150mL



11.8% More milk¹
105° angled breast shield
provides comfort on
the breast



Comfortable soft rim



Closed system
prevents breast milk from entering
the tubing and motor



Oval for better fit



Big on Performance. Loved by Moms.

10
OUT OF
10

100% of moms are **highly satisfied** with the **milk output***

9
OUT OF
10

98% of moms said Medela Freestyle™ Mini Hands-free is **ultra small but powerful***

10
OUT OF
10

100% of moms found expressing breast milk with the Medela Freestyle™ Mini Hands-free breast pump **very comfortable***



What's included

- x1 Freestyle™ Mini Breast Pump
- x2 Hands-free Collection Cups with 24mm breast shield and O-ring
- x2 3 sizes of inserts:
15mm / 18mm / 21mm
- x1 Adjustable lanyard with clip-on feature
- x1 Double tubing
- x1 Pouch
- x1 USB-C cable
- x4 Easy Pour Breast Milk Storage Bags

Item Number	Product Name	UPC	Package Information			Shipper Information		
			Unit Count	Dimensions (in.)	Weight (lbs.)	Qty/Case	Dimensions (in.)	Weight (lbs.)
101048395	Freestyle™ Mini Hands-free	020451904870	1	4.92 x 8.66 x 5.62	2.2	1	5.51 x 9.01 x 6.10	2.2



Medela wordmark and logo are registered in the U.S. Patent and Trademark Office and the Canadian Intellectual Property Office. Freestyle is a trademark of Medela. 0226 © 2026 Medela

*Based on the 10 days usage in the in-home study in Mexico with 50 mothers, Sept 2025.
**Applies when using without the insert.
1 Sakalidis VS et al. Acta Obstet Gynecol Scand. 2020; 99(11):1561-1567.