

Improve your quality of life with Portable Oxygen Conctrators

A portable oxygen concentrator (**POC**) is designed to use the air around you and turn it into medical oxygen prescribed by your doctor.

Lightweight & hands-free with various carry case options, **POC's** are designed to move with you!

Forget the days of lugging around heavy oxygen cylinders and breathe easy, a **POC** is made for your on-the-go lifestyle.



Scan the QR code or visit carelincmcd.com/locations to find the closest CareLinc location to learn more!



Every OxyGo NEXT comes with:

- Single battery
- Nasal cannula
- AC & DC chargers
- Carry bag with padded adjustable strap



Every OxyGo FIT comes with:

- A durable, black protective bag
- Carry strap
- AC & DC chargers
- Easy to follow instructions



OxyGoNEXT™

OxyGoNEXT™

OxyGoFIT™

OxyGoFIT™

	Single Battery	Double Battery	Single Battery	Double Battery
Weight	4.8 lbs	5.8 lbs	2.8 lbs	3.29 lbs
Size	7.18" x 3.27" x 8.14"	7.18" x 3.27" x 9.02"	5.91" x 2.68" x 7.2"	5.91" x 2.68" x 7.79"
Battery Duration	Up to 6.5 hours	Up to 12.75 hours	Up to 2 hours	Up to 5 hours
Recharge Time	Up to 3 hours	Up to 4 hours	Up to 2.7 hours	Up to 5 hours
Flow Settings	1 (210 ml/min) 2 (420 ml/min) 3 (630 ml/min) 4 (840 ml/min) 5 (1050 ml/min) 6 (1260 ml/min)		1 (210 ml/min) 2 (420 ml/min) 3 (630 ml/min)	
Power	AC & DC power supply		AC & DC power supply	
Warranty	5 year limited		3 year limited	
Sound Level (at flow 2)	38 decibels		40 decibels	