



Medical Equipment & Supply

www.carelinmed.com

Wheelchair Sizing Chart

by **Height** and **Weight**

***Note: This chart is approximate, please measure for exact fit.**

	< 5'0"	5'0"	5'2"	5'4"	5'6"	5'8"	5'10"	6'0"	6'2"	6'4"	6'6"	> 6'6"
< 100lbs	16"	16"	16"	16"	16"	16"	16"	16"	16"	16"	16"	16"
100lbs	16"	16"	16"	16"	16"/18"	16"	16"	16"	16"	16"	16"	16"
120lbs	18"	18"	18"	18"	16"/18"	16"/18"	16"	16"	16"	16"	16"	16"
140lbs	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"	18"	16"	16"
160lbs	20"	20"	18"/20"	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"
180lbs	20"	20"	20"	20"	20"	18"/20"	18/20"	18"	18"	18"	18"	18"
200lbs	22"	22"	20"	20"	20"	20"	20"	18"/20"	18"/20"	18"/20"	18"	18"
220lbs	22"	22"	22"	20"	20"	20"	20"	20"	20"	20"	20"	20"
240lbs	24"	24"	22"/24"	22"/24"	20"	20"	22"	22"	22"	22"	22"	22"
260lbs	24"	24"	24"	24"	22"/24"	22"/24"	22"	22"	22"	22"	22"	22"
280lbs	CONTACT CARELINC			24"	24"	24"	22"/24"	22"	22"	22"	22"	22"
300lbs						24"	24"	24"	24"	22"	22"	22"
320lbs							24"	24"	24"	24"	24"	24"
340lbs								24"	24"	24"	24"	24"
								24"	24"	24"	24"	24"

Seat Length:
16"
18"
20"

Any patient under 5'2" - Hemi height, patients 5'4" should use 2" cushion. Patients in-between sizes, go down for men and up for women. If you see "Contact CareLinc" in your size range, a standard chair may not be the right fit. Call us or visit carelinmed.com so our team can help determine the right chair for the patient.