



Power Mobility Battery Tips

Stay Charged & Ready to Roll!



A Little Maintenance Goes a Long Way

Just like your car needs an oil change, your power chair needs a yearly checkup. We recommend annual maintenance—and your batteries might need to be swapped out once a year too. Even with perfect care, no battery lasts forever! Follow these battery tips to make sure your equipment is always ready to go when you are.

Charging Guidelines:

- Daily driver? Charge it every night for 8–14 hours.
- Occasional outings? Charge it 2–3 times a week for 8–14 hours each time.
- Storing for more than 2 weeks? Give it a full charge, then disconnect the batteries.
- At the very least, plug in and charge your chair once a month for 8–14 hours.
- Keep your equipment performing at its best with an annual checkup, covered by insurance.

Things to Avoid:

- Don't overdo it! Charging for more than 24 hours straight can wear things out.
- Don't run it dry! Avoid draining your battery all the way down to the last red light.
- No car chargers allowed! Only use the charger that came with your equipment—no wet-cell or car battery chargers.
- Keep the battery away from areas with excessive moisture or freezing temperatures.

Need Help or Have Questions?

Phone: 616.249.2273, option 5
Email: rehabservice@carelinmed.com
www.carelinmed.com



Ready to make your home safer and more accessible?
Scan the QR code or visit carelinmed.click/home-access
to explore solutions designed around your needs.