







All our breast pumps are
closed-systems & BPA free!

Start your order at
www.carelincm.com/breast-pumps

Insurance Billable Breast Pumps Comparison Chart

							
WITH INSURANCE		\$0	\$0	\$0	\$0	UPGRADE \$135	UPGRADE \$200
BRAND		Medela	Medela	Motif	Spectra	Spectra	Medela
MODEL		Pump in Style Pro	Swing Hands-Free	Luna with Battery	S2 Plus	S1 Plus	Freestyle Hands-Free
WEIGHT		1.8 lbs	< 1 lb	1.9 lbs	3 lbs	3 lbs	< 1 lb
SUCTION STRENGTH		270 mmHg	245 mmHg	280 mmHg	270 mmHg	270 mmHg	245 mmHg
WARRANTY		1 yr. Motor 90 day parts	1 yr. Motor 90 day parts	2 yr. Motor 90 day parts	2 yr. Motor 90 day parts	2 yr. Motor 90 day parts	1 yr. Motor 90 day parts
BATTERY POWER		AA-8 (not included)	✓	✓	✗	✓	✓
AUTO-SHUTOFF		✗	✓	✓	✓	✓	✓
DIGITAL DISPLAY		✗	✓	✓	✓	✓	✓
BACKLIT DISPLAY		✓	✓	✓	✗	✓	✓
LAST SETTING MEMORY		✗	✓	✓	✓	✓	✓
EXTRAS		Pumping similar to Symphony®, 16 Levels, Mode Indicator, Level 3 start	Hands-Free, Rechargeable Battery	Small but powerful, Rechargeable Battery	Nightlight, LED Timer	Same as S2 Plus, w/ Rechargeable Battery, Cooler, & Tote Bag	Hands-Free, Rechargeable Battery

* Offered exclusively online at carelincm.com in our SHOP & ships direct to home. Cannot be purchased at your local CareLinc Medical store.

SS11CE00308

CARELINC