

Insurance-Billable Breast Pumps

All breast pump options from CareLinc are double electric & BPA free to protect the health of you and your child!



Specs	Advanced Pump	Pump-in-Style	S2Plus
Description	Evenflo's Advanced Pump features the PerfectlyPure™ design to prevent backflow & cleans up easily. AdvancedControl™ technology with 32 speed & suction settings.	Medela's Pump-in-Style Advanced kit delivers effective pumping in a modular solution, offering 2-Phase Expression technology.	The Spectra S2Plus is a pump designed for comfortability and includes features such as massage mode & back-flow protection, LCD display with timer & nightlight.
What's Included	2 each - 24.5mm/28mm Soft Flange Inserts 2 - 30.5mm Flange Kits 2- 5oz Storage Bottles 2 - Extra Valves & Membranes 1 - AC Adapter Documents: Instructions for use & to access online breastfeeding & breast pump instruct programs, breastfeeding guide card	2 each - 24mm/27mm Breast Shields 2 - Breast Shield Connectors 2 - Valves & Membranes 2 - 5 oz Storage Bottles 1 - Complete Tubing 1 - AC Adapter Documents: Instructions for Use	2 - each 24mm/28mm Flanges 2 - Tubing 2 - Backflow Protectors 2 - Duckbill Valves 2 - 5 oz Storage Bottles 1 - AC Adapter Documents: Product manual
Carrying Bag	None Available Through CareLinc	*See Upgrade	*See Upgrade
Settings	Independent Suction/Speed Adjustments	Dial For Simultaneous Adjustments	Independent Cycle/Vacuum Adjustments
Closed System	Yes	No	Yes
mmHg	Up to 250 mmHg	Up to 270 mmHg	Up to 280 mmHg
Battery	Sold Separately	Sold Separately	*See Upgrade
Weight / Dimensions	≈ 3 lbs. / 12" x 5.6" x 9.69"	≈ 5 lbs. / 9.25" x 7.5" x 8.5"	≈ 2.5 lbs. / 7.9" x 7.9" x 6.7"
Power	Dual Options	Dual Options	AC Adapter Only
Warranty	1 Year on Motor, 90 days for parts	1 Year on Motor, 90 days for parts	2 Years on Motor, 90 days for parts
*Upgrade	Not Available Through CareLinc	On-The-Go Tote (\$115)	Spectra S1Plus / Cooler & Tote (\$115)

The type of breast pump available to you may depend on your insurance coverage. We are contracted with most insurance plans. Additional breastfeeding accessories are available for purchase in-store & at www.carelincmed.com.

* Upgrade available with insurance coverage

\$115



medela
On-The-Go Tote



speCtra®
S1Plus

Start your order at www.carelincmed.com